

## Reptin / RUVBL2 Rabbit mAb

Catalog No: #56203

Package Size: #56203-1 50ul #56203-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Reptin / RUVBL2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse
Specificity	Reptin / RUVBL2 Antibody detects endogenous levels of total Reptin / RUVBL2
Immunogen Description	A synthesized peptide derived from human Reptin / RUVBL2
Other Names	ECP51; Reptin 52; REPTIN; RUVB2; RUVBL2; RVB2; TAP54-beta; TIH2; TIP48; TIP49b; zreptin;
Accession No.	Uniprot:Q9Y230
Uniprot	Q9Y230
Calculated MW	48kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of Reptin / RUVBL2 expression in (1) HeLa cell lysate; (2) NIH/3T3 cell lysate.

## Product Description

Possesses single-stranded DNA-stimulated ATPase and ATP-dependent DNA helicase (5' to 3') activity. Component of the NuA4 histone acetyltransferase complex which is involved in transcriptional activation of select genes principally by acetylation of nucleosomal histones H4 and H2A.

---

Note: This product is for in vitro research use only