

ORM1 Rabbit mAb

Catalog No: #56212

Package Size: #56212-1 50ul #56212-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ORM1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	ORM1 Antibody detects endogenous levels of total ORM1
Immunogen Description	A synthesized peptide derived from human ORM1
Other Names	AGP 1; AGP A; AGP1; Alpha 1 acid glycoprotein; alpha-1-AGP; HEL S 153w; ORM; Orosomucoid 1;
Accession No.	Uniprot:P02763
Uniprot	P02763
Calculated MW	47kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of ORM1 expression in human testis lysate.

Product Description

Appears to function in modulating the activity of the immune system during the acute-phase reaction.

Note: This product is for in vitro research use only