

INDO Rabbit mAb

Catalog No: #56213

Package Size: #56213-1 50ul #56213-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	INDO Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human
Specificity	INDO Antibody detects endogenous levels of total INDO
Immunogen Description	A synthesized peptide derived from human INDO
Other Names	3-dioxygenase; IDO 1; IDO; IDO1; INDO; indolamine 2,3 dioxygenase; Indole 2 3 dioxygenase;
Accession No.	Uniprot:P14902
Uniprot	P14902
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000 ICC/IF:1:50~1:200

Images



Western blot analysis of INDO expression in HeLa cell lysate treated with IFN gamma.

Product Description

Catalyzes the cleavage of the pyrrol ring of tryptophan and incorporates both atoms of a molecule of oxygen.

Note: This product is for in vitro research use only