

ADORA1 Rabbit mAb

Catalog No: #56216

Package Size: #56216-1 50ul #56216-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ADORA1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	ADORA1 Antibody detects endogenous levels of total ADORA1
Immunogen Description	A synthesized peptide derived from human ADORA1
Other Names	A1AR; A1R; AA1R; Adenosine receptor A1; ADORA 1; ADORA1; RDC 7; Ri;
Accession No.	Uniprot:P30542
Uniprot	P30542
Calculated MW	37kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50FC:1:100

Images



Western blot analysis of ADORA1 expression in A549 cell lysate.

Product Description

Receptor for adenosine. The activity of this receptor is mediated by G proteins which inhibit adenylyl cyclase.

Note: This product is for in vitro research use only