

## Grp75 Rabbit mAb

Catalog No: #56225

Package Size: #56225-1 50ul #56225-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Grp75 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Grp75 Antibody detects endogenous levels of total Grp75
Immunogen Description	A synthesized peptide derived from human Grp75
Other Names	GRP-75; Heat shock 70 kDa protein 9; Mortalin; MOT; Peptide-binding protein 74; PBP74; HSPA9; GRP75; HSPA9B; mt-HSP70;
Accession No.	Uniprot:P38646
Uniprot	P38646
Calculated MW	74kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.


## Application Details

WB:1:500~1:2000

## Images



Western blot analysis of Grp75 expression in (1) HeLa cell lysate; (2) Mouse liver lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

## Product Description

Implicated in the control of cell proliferation and cellular aging. May also act as a chaperone.

Note: This product is for in vitro research use only