

PDE4D Rabbit mAb

Catalog No: #56238

Package Size: #56238-1 50ul #56238-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PDE4D Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	PDE4D Antibody detects endogenous levels of total PDE4D
Immunogen Description	A synthesized peptide derived from human PDE4D
Other Names	cAMP-specific 3',5'-cyclic phosphodiesterase 4D; DPDE3; PDE43; PDE4D; DPDE3;
Accession No.	Uniprot:Q08499
Uniprot	Q08499
Calculated MW	91kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of PDE4D expression in human fetal heart lysate.

Product Description

Hydrolyzes the second messenger cAMP, which is a key regulator of many important physiological processes.

Note: This product is for in vitro research use only