

COUP TF1 Rabbit mAb

Catalog No: #56245

Package Size: #56245-1 50ul #56245-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | COUP TF1 Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB IHC ICC/IF FC |
| Species Reactivity | Human Mouse Rat |
| Specificity | COUP TF1 Antibody detects endogenous levels of total COUP TF1 |
| Immunogen Description | A synthesized peptide derived from human COUP TF1 |
| Other Names | COUP-TF1; EAR3; ERBAL3; NR2F1; NR2F2; Nuclear receptor subfamily 2 group F member 1; SVP44; TCFCOUP1; TFCOUP1; |
| Accession No. | Uniprot:P10589 |
| Uniprot | P10589 |
| Calculated MW | 46kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:500~1:1000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of COUP TF1 expression in HEK293 cell lysate.

Product Description

Coup (chicken ovalbumin upstream promoter) transcription factor binds to the ovalbumin promoter and, in conjunction with another protein (S300-II) stimulates initiation of transcription. Binds to both direct repeats and palindromes of the 5'-AGGTCA-3' motif. Represses transcriptional activity of LHCG.

Note: This product is for in vitro research use only