

# Trefoil Factor 3 Rabbit mAb

Catalog No: #56250

Package Size: #56250-1 50ul #56250-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Trefoil Factor 3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	Trefoil Factor 3 Antibody detects endogenous levels of total Trefoil Factor 3
Immunogen Description	A synthesized peptide derived from human Trefoil Factor 3
Other Names	hITF; hP1.B; Intestinal trefoil factor; ITF; mITF; P1B; Polypeptide P1.B; TFF3; TFI; TREFOIL; Trefoil factor 3;
Accession No.	Uniprot:Q07654
Uniprot	Q07654
Calculated MW	9kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:1000IHC:1:100~1:500IP:1:30

## Images



Western blot analysis of Trefoil Factor 3 expression in human small intestine lysate.

## Product Description

Involved in the maintenance and repair of the intestinal mucosa. Promotes the mobility of epithelial cells in healing processes (motogen).

---

Note: This product is for in vitro research use only