

## CD82 Rabbit mAb

Catalog No: #56256

Package Size: #56256-1 50ul #56256-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	CD82 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human Rat
Specificity	CD82 Antibody detects endogenous levels of total CD82
Immunogen Description	A synthesized peptide derived from human CD82
Other Names	C33; CD82; IA4; Inducible membrane protein; KAI1; Kangai1; SAR2; ST6; Tetraspanin 27; Tspan27;
Accession No.	Uniprot:P27701
Uniprot	P27701
Calculated MW	30-60kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.


## Application Details

WB: 1:500~1:2000

IHC: 1:50~1:200


FC: 1:50

## Images



Western blot analysis of CD82 expression in Jurkat cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



## Product Description

Associates with CD4 or CD8 and delivers costimulatory signals for the TCR/CD3 pathway.

Note: This product is for in vitro research use only