

TRADD Rabbit mAb

Catalog No: #56274

Package Size: #56274-1 50ul #56274-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TRADD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC IP
Species Reactivity	Human Mouse Rat
Specificity	TRADD Antibody detects endogenous levels of total TRADD
Immunogen Description	A synthesized peptide derived from human TRADD
Other Names	AA930854; TNFR1 associated DEATH domain protein; tradd;
Accession No.	Uniprot:Q15628
Uniprot	Q15628
Calculated MW	34kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:1000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:10

Images



Western blot analysis of TRADD expression in Hela cell lysate.

Product Description

Adapter molecule for TNFRSF1A/TNFR1 that specifically associates with the cytoplasmic domain of activated TNFRSF1A/TNFR1 mediating its interaction with FADD. Overexpression of TRADD leads to two major TNF-induced responses, apoptosis and activation of NF-kappa-B.

Note: This product is for in vitro research use only