

CD33 Rabbit mAb

Catalog No: #56297

Package Size: #56297-1 50ul #56297-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD33 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	CD33 Antibody detects endogenous levels of total CD33
Immunogen Description	A synthesized peptide derived from human CD33
Other Names	CD33; gp67; My9; p67; SIGLEC3;
Accession No.	Uniprot:P20138
Uniprot	P20138
Calculated MW	65-75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:50

Images



Western blot analysis of CD33 expression in THP1 cell lysate.

Product Description

Putative adhesion molecule of myelomonocytic-derived cells that mediates sialic-acid dependent binding to cells. Preferentially binds to alpha-2,6-linked sialic acid.

Note: This product is for in vitro research use only