

alpha Tubulin (yeast) Rabbit mAb

Catalog No: #56316

Package Size: #56316-1 50ul #56316-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	alpha Tubulin (yeast) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Saccharomyces cerevisiae yeast
Specificity	alpha Tubulin yeast Antibody detects endogenous levels of total alpha Tubulin yeast
Immunogen Description	A synthesized peptide derived from yeast alpha Tubulin
Other Names	Alpha-tubulin 1;H2 ALPHA;LIS3;TUBA1;TUBA1A;Tuba4a;Tubulin alpha-4A chain
Accession No.	Uniprot:P09733
Uniprot	P09733
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:3000~1:10000

Images



Western blot analysis of alpha Tubulin expression in *Saccharomyces cerevisiae* (yeast) cell lysate.

Product Description

Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.

Note: This product is for in vitro research use only