

## FHL1 Rabbit mAb

Catalog No: #56326

Package Size: #56326-1 50ul #56326-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	FHL1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	FHL1 Antibody detects endogenous levels of total FHL1
Immunogen Description	A synthesized peptide derived from human FHL1
Other Names	FHL1; FHL1A; FHL1B; FLH1A; RBMX1A; RBMX1B; SLIM1; SLIMMER; XMPMA;
Accession No.	Uniprot:Q13642
Uniprot	Q13642
Calculated MW	30kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

## Images



Western blot analysis of FHL1 expression in Human fetal muscle lysate.

## Product Description

May have an involvement in muscle development or hypertrophy.

---

Note: This product is for in vitro research use only