

Cytochrome P450 1A2 Rabbit mAb

Catalog No: #56329

Package Size: #56329-1 50ul #56329-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Cytochrome P450 1A2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	Cytochrome P450 1A2 Antibody detects endogenous levels of total Cytochrome P450 1A2
Immunogen Description	A synthesized peptide derived from human Cytochrome P450 1A2
Other Names	CP12; CYP1A2; CYP1A2; P3 450; P450 4; P450 P3;
Accession No.	Uniprot:P05177
Uniprot	P05177
Calculated MW	58kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200FC:1:500

Images



Western blot analysis of Cytochrome P450 1A2 expression in Caco2 cell lysate.

Product Description

Cytochromes P450 are a group of heme-thiolate monooxygenases. In liver microsomes, this enzyme is involved in an NADPH-dependent electron transport pathway.

Note: This product is for in vitro research use only