

## Thyroglobulin Rabbit mAb

Catalog No: #56331

Package Size: #56331-1 50ul #56331-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Thyroglobulin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	Thyroglobulin Antibody detects endogenous levels of total Thyroglobulin
Immunogen Description	A synthesized peptide derived from human Thyroglobulin
Other Names	AITD3; hTG; TDH3; Tg; TGN; Thyroglobulin;
Accession No.	Uniprot:P01266
Uniprot	P01266
Calculated MW	305kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:1000~1:5000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

## Images



Western blot analysis of Thyroglobulin expression in Mouse thyroid lysate.

## Product Description

Precursor of the iodinated thyroid hormones thyroxine (T4) and triiodothyronine (T3).

---

Note: This product is for in vitro research use only