

beta Defensin 1 Rabbit mAb

Catalog No: #56333

Package Size: #56333-1 50ul #56333-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	beta Defensin 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	beta Defensin 1 Antibody detects endogenous levels of total beta Defensin 1
Immunogen Description	A synthesized peptide derived from human beta Defensin 1
Other Names	BD1; beta 1; DEFB1; DEFB101; Defensin; Defensin beta 1; HBD1;
Accession No.	Uniprot:P60022
Uniprot	P60022
Calculated MW	7kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of beta Defensin 1 expression in Human recombinant beta Defensin 1 lysate.

Product Description

Has bactericidal activity.

Note: This product is for in vitro research use only