

Mammaglobin A Rabbit mAb

Catalog No: #56348

Package Size: #56348-1 50ul #56348-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Mammaglobin A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Rat
Specificity	Mammaglobin A Antibody detects endogenous levels of total Mammaglobin A
Immunogen Description	A synthesized peptide derived from human Mammaglobin A
Other Names	Mammaglobin 1; Mammaglobin-A; MGB1; Scgb2a2; SG2A2; UGB2;
Accession No.	Uniprot:Q13296
Uniprot	Q13296
Calculated MW	18kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200

Images



Western blot analysis of Mammaglobin A expression in Human breast cancer lysate.

Product Description

Belongs to the secretoglobulin family. Lipophilin subfamily.

Note: This product is for in vitro research use only