

## DMRT1 Rabbit mAb

Catalog No: #56354

Package Size: #56354-1 50ul #56354-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	DMRT1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	DMRT1 Antibody detects endogenous levels of total DMRT1
Immunogen Description	A synthesized peptide derived from human DMRT1
Other Names	CT154; dmr1; DMT1;
Accession No.	Uniprot:Q9Y5R6
Uniprot	Q9Y5R6
Calculated MW	34kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IP:1:80

## Images



Western blot analysis of DMRT1 expression in HL-60 cell lysate.

## Product Description

May be required for testis development.

---

Note: This product is for in vitro research use only