

## MSMB Rabbit mAb

Catalog No: #56355

Package Size: #56355-1 50ul #56355-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	MSMB Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	Prostate Secretory Protein/PSP Antibody detects endogenous levels of total Prostate Secretory Protein/PSP
Immunogen Description	A synthesized peptide derived from human Prostate Secretory Protein/PSP
Other Names	MSP; PSP; IGBF; MSPB; PN44; PRPS; HPC13; PSP57; PSP94; PSP-94;
Accession No.	Uniprot:P08118
Uniprot	P08118
Calculated MW	15kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:1000~1:5000IHC:1:100~1:500IP:1:50

## Images



Western blot analysis of Prostate Secretory Protein/PSP expression in Human prostate cancer lysate.

## Product Description

Specific receptors for this protein are found on spermatozoa and in the prostate.

---

Note: This product is for in vitro research use only