

CD79B Rabbit mAb

Catalog No: #56362

Package Size: #56362-1 50ul #56362-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CD79B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	CD79b Antibody detects endogenous levels of total CD79b
Immunogen Description	A synthesized peptide derived from human CD79b
Other Names	CD79B; AGM6; B29; IGB;
Accession No.	Uniprot:P40259
Uniprot	P40259
Calculated MW	30-50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50FC:1:50

Images



Western blot analysis of CD79b expression in Raji cell lysate.

Product Description

Required in cooperation with CD79A for initiation of the signal transduction cascade activated by the B-cell antigen receptor complex (BCR) which leads to internalization of the complex, trafficking to late endosomes and antigen presentation.

Note: This product is for in vitro research use only