

## PMP70 Rabbit mAb

Catalog No: #56364

Package Size: #56364-1 50ul #56364-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PMP70 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	PMP70 Antibody detects endogenous levels of total PMP70
Immunogen Description	A synthesized peptide derived from human PMP70
Other Names	ABCD3; ABC43; PMP70; PXMP1; ZWS2;
Accession No.	Uniprot:P28288
Uniprot	P28288
Calculated MW	70kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200FC:1:100

## Images



Western blot analysis of PMP70 expression in HepG2 cell lysate.

## Product Description

Probable transporter. The nucleotide-binding fold acts as an ATP-binding subunit with ATPase activity.

---

Note: This product is for in vitro research use only