

## RANK Rabbit mAb

Catalog No: #56372

Package Size: #56372-1 50ul #56372-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	RANK Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	RANK Antibody detects endogenous levels of total RANK
Immunogen Description	A synthesized peptide derived from human RANK
Other Names	CD265; FEO; LOH18CR1; mRANK; ODFR; OFE; OPTB7; OST5; PDB2; RANK; TNFRSF11A; TRANCER;
Accession No.	Uniprot:Q9Y6Q6
Uniprot	Q9Y6Q6
Calculated MW	80kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of RANK expression in Human fetal liver lysate.

## Product Description

Receptor for TNFSF11/RANKL/TRANCE/OPGL; essential for RANKL-mediated osteoclastogenesis. Involved in the regulation of interactions between T-cells and dendritic cells.

---

Note: This product is for in vitro research use only