

OLFM4 Rabbit mAb

Catalog No: #56375

Package Size: #56375-1 50ul #56375-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	OLFM4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	OLFM4 Antibody detects endogenous levels of total OLFM4
Immunogen Description	A synthesized peptide derived from human OLFM4
Other Names	GC1; GW112; hGC 1; hOLFID; olfactomedin 4; OlfD; OLFM4; OLM4;
Accession No.	Uniprot:Q6UX06
Uniprot	Q6UX06
Calculated MW	54kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:40

Images



Western blot analysis of OLFM4 expression in LnCaP cell lysate.

Product Description

An antiapoptotic factor that promotes tumor growth. Promotes proliferation of pancreatic cancer cells by favoring the transition from the S to G2/M phase. Facilitates cell adhesion.

Note: This product is for in vitro research use only