

## Alkaline phosphatase Rabbit mAb

Catalog No: #56377

Package Size: #56377-1 50ul #56377-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	Alkaline phosphatase Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	intestinal alkaline phosphatase Antibody detects endogenous levels of total intestinal alkaline phosphatase
Immunogen Description	A synthesized peptide derived from human intestinal alkaline phosphatase
Other Names	Alpi; Glycerophosphatase; IAP;
Accession No.	Uniprot:P09923
Uniprot	P09923
Calculated MW	68kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000


IHC: 1:50~1:200

## Images



Western blot analysis of intestinal alkaline phosphatase expression in HeLa cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



## Product Description

In most mammals there are four different isozymes: placental, placental-like, intestinal and tissue non-specific (liver/bone/kidney).

Note: This product is for in vitro research use only