

CD53 Rabbit mAb

Catalog No: #56387

Package Size: #56387-1 50ul #56387-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD53 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	CD53 Antibody detects endogenous levels of total CD53
Immunogen Description	A synthesized peptide derived from human CD53
Other Names	CD53; CD53 antigen; MOX44; OX44; RATOX44; Tetraspanin25; Transmembrane glycoprotein; TSPAN25; PRO1340;
Accession No.	Uniprot:P19397
Uniprot	P19397
Calculated MW	24kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:100~1:500FC:1:50

Images



Western blot analysis of CD53 expression in Jurkat cell lysate.

Product Description

May be involved in growth regulation in hematopoietic cells.

Note: This product is for in vitro research use only