

Myogenin Rabbit mAb

Catalog No: #56388

Package Size: #56388-1 50ul #56388-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | Myogenin Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB IHC FC |
| Species Reactivity | Human Mouse Rat |
| Specificity | Myogenin Antibody detects endogenous levels of total Myogenin |
| Immunogen Description | A synthesized peptide derived from human Myogenin |
| Other Names | bHLHc3; cb553; MYF4; MYOG; Myogenic factor 4; |
| Accession No. | Uniprot:P15173 |
| Uniprot | P15173 |
| Calculated MW | 36kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:500~1:2000IHC:1:50~1:200FC:1:40

Images



Western blot analysis of Myogenin expression in C2C12 cell lysate.

Product Description

Involved in muscle differentiation (myogenic factor). Induces fibroblasts to differentiate into myoblasts. Probable sequence specific DNA-binding protein.

Note: This product is for in vitro research use only