

TPO Rabbit mAb

Catalog No: #56389

Package Size: #56389-1 50ul #56389-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name | TPO Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB IHC IP |
| Species Reactivity | Human |
| Specificity | TPO Antibody detects endogenous levels of total TPO |
| Immunogen Description | A synthesized peptide derived from human TPO |
| Other Names | MSA; TDH2A; TPO; TPX; |
| Accession No. | Uniprot:P07202 |
| Uniprot | P07202 |
| Calculated MW | 103kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:40

Images



Western blot analysis of TPO expression in Human thyroid lysate.

Product Description

Iodination and coupling of the hormonogenic tyrosines in thyroglobulin to yield the thyroid hormones T(3) and T(4).

Note: This product is for in vitro research use only