

SCGN Rabbit mAb

Catalog No: #56393

Package Size: #56393-1 50ul #56393-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SCGN Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	SCGN Antibody detects endogenous levels of total SCGN
Immunogen Description	A synthesized peptide derived from human SCGN
Other Names	CALBL; Scgn; SECRET; Secretagogin; SEGN; Setagin;
Accession No.	Uniprot:O76038
Uniprot	O76038
Calculated MW	32kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of SCGN expression in human fetal brain lysate.

Product Description

Expressed at high levels in the pancreatic islets of Langerhans and to a much lesser extent in the gastrointestinal tract (stomach, small intestine and colon), the adrenal medulla and cortex and the thyroid C-cells.

Note: This product is for in vitro research use only