

## Thyroglobulin Rabbit mAb

Catalog No: #56397

Package Size: #56397-1 50ul #56397-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

|                       |  |
|-----------------------|--|
| Product Name          | Thyroglobulin Rabbit mAb   |
| Host Species          | Rabbit   |
| Clonality             | Monoclonal   |
| Isotype               | Rabbit IgG   |
| Purification          | Affinity-chromatography  |
| Applications          | WB IHC IP  |
| Species Reactivity    | Human  |
| Specificity           | Thyroglobulin Antibody detects endogenous levels of total Thyroglobulin                            |
| Immunogen Description | A synthesized peptide derived from human Thyroglobulin   |
| Other Names           | AITD3; hTG; TDH3; Tg; TGN; Thyroglobulin;  |
| Accession No.         | Uniprot:P01266   |
| Uniprot               | P01266   |
| Calculated MW         | 305kDa   |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle. |

## Application Details

WB:1:1000~1:5000IHC:1:50~1:200IP:1:50

## Images



Western blot analysis of Thyroglobulin expression in Human thyroid lysate.

## Product Description

Precursor of the iodinated thyroid hormones thyroxine (T4) and triiodothyronine (T3).

---

Note: This product is for in vitro research use only