

TBX21 Rabbit mAb

Catalog No: #56399

Package Size: #56399-1 50ul #56399-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TBX21 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IHC FC
Species Reactivity	Human
Specificity	TBX21 Antibody detects endogenous levels of total TBX21
Immunogen Description	A synthesized peptide derived from human TBX21
Other Names	T bet; T box 21; T PET; TBET; TBLYM; Tbx21; TPET; Transcription factor TBLYM;
Accession No.	Uniprot:Q9UL17
Uniprot	Q9UL17
Calculated MW	kDa
Storage	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

IHC:1:50~1:200FC:1:50

Images



Immunohistochemical analysis of paraffin-embedded human spleen, using TBX21 Antibody.

Product Description

Transcription factor that controls the expression of the TH1 cytokine, interferon-gamma. Initiates TH1 lineage development from naive TH precursor cells both by activating TH1 genetic programs and by repressing the opposing TH2 programs.

Note: This product is for in vitro research use only