

PBRM1 Rabbit mAb

Catalog No: #56406

Package Size: #56406-1 50ul #56406-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PBRM1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	PBRM1 Antibody detects endogenous levels of total PBRM1
Immunogen Description	A synthesized peptide derived from human PBRM1
Other Names	BAF180; hPB1; PB1; Pbrm1; polybromo 1; Polybromo-1D;
Accession No.	Uniprot:Q86U86
Uniprot	Q86U86
Calculated MW	193kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:100~1:500ICC/IF:1:50~1:100FC:1:50

Images



Western blot analysis of PBRM1 expression in HeLa cell lysate.

Product Description

Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-nucleosome topology). Act as a negative regulator of cell proliferation.

Note: This product is for in vitro research use only