

TRBP Rabbit mAb

Catalog No: #56410

Package Size: #56410-1 50ul #56410-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TRBP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	TRBP Antibody detects endogenous levels of total TRBP
Immunogen Description	A synthesized peptide derived from human TRBP
Other Names	LOQS; Prbp; tarbp2; TRBP; TRBP1; TRBP2;
Accession No.	Uniprot:Q15633
Uniprot	Q15633
Calculated MW	39kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:30FC:1:50

Images



Western blot analysis of TRBP expression in Jurkat cell lysate.

Product Description

Required for formation of the RNA induced silencing complex (RISC). Component of the RISC loading complex (RLC), also known as the micro-RNA (miRNA) loading complex (miRLC), which is composed of DICER1, EIF2C2/AGO2 and TARBP2.

Note: This product is for in vitro research use only