

HES5 Rabbit mAb

Catalog No: #56412

Package Size: #56412-1 50ul #56412-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HES5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	HES5 Antibody detects endogenous levels of total HES5
Immunogen Description	A synthesized peptide derived from human HES5
Other Names	HES 5; Hes5;
Accession No.	Uniprot:Q5TA89
Uniprot	Q5TA89
Calculated MW	18kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of HES5 expression in (1) 293 cell lysate; (2) Mouse spleen lysate; (3) C6 cell lysate.

Product Description

Transcriptional repressor of genes that require a bHLH protein for their transcription.

Note: This product is for in vitro research use only