

## ALDH2 Rabbit mAb

Catalog No: #56421

Package Size: #56421-1 50ul #56421-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ALDH2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	ALDH2 Antibody detects endogenous levels of total ALDH2
Immunogen Description	A synthesized peptide derived from human ALDH2
Other Names	ALDH 2; ALDH E2; Aldh2; ALDH1; ALDM; MS767;
Accession No.	Uniprot:P05091
Uniprot	P05091
Calculated MW	56kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50

## Images



Western blot analysis of ALDH2 expression in HepG2 cell lysate.

## Product Description

Alcohol metabolism; ethanol degradation; acetate from ethanol: step 2/2.

---

Note: This product is for in vitro research use only