

CD134 Rabbit mAb

Catalog No: #56422

Package Size: #56422-1 50ul #56422-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CD134 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	CD134 Antibody detects endogenous levels of total CD134
Immunogen Description	A synthesized peptide derived from human CD134
Other Names	ACT35; CD134; IMD16; OX40; OX40L receptor; TNFRSF4; Txgp1; Txgp1l;
Accession No.	Uniprot:P43489
Uniprot	P43489
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:1000FC:1:50

Images



Western blot analysis of CD134 expression in Human fetal spleen lysate.

Product Description

Receptor for TNFSF4/OX40L/GP34.

Note: This product is for in vitro research use only