

## E2F6 Rabbit mAb

Catalog No: #56426

Package Size: #56426-1 50ul #56426-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	E2F6 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	E2F6 Antibody detects endogenous levels of total E2F6
Immunogen Description	A synthesized peptide derived from human E2F6
Other Names	E2F 6; E2F-6; E2F6; EMA;
Accession No.	Uniprot:O75461
Uniprot	O75461
Calculated MW	32kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000 ICC/IF:1:50~1:200 FC:1:50

## Images



Western blot analysis of E2F6 expression in K562 cell lysate.

## Product Description

Inhibitor of E2F-dependent transcription. Binds DNA cooperatively with DP proteins through the E2 recognition site, 5'-TTTC[CG]CGC-3'. Has a preference for the 5'-TTTCCCGC-3' E2F recognition site. E2F-6 lacks the transcriptional activation and pocket protein binding domains.

---

Note: This product is for in vitro research use only