

CD89 Rabbit mAb

Catalog No: #56433

Package Size: #56433-1 50ul #56433-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CD89 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	CD89 Antibody detects endogenous levels of total CD89
Immunogen Description	A synthesized peptide derived from human CD89
Other Names	CD89; Fc alpha receptor; FCAR; IgA Fc receptor;
Accession No.	Uniprot:P24071
Uniprot	P24071
Calculated MW	70kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:100~1:500ICC/IF:1:100~1:500IP:1:50FC:1:50

Images



Western blot analysis of CD89 expression in HL-60 cell lysate.

Product Description

Binds to the Fc region of immunoglobulins alpha. Mediates several functions including cytokine production.

Note: This product is for in vitro research use only