

CD66b Rabbit mAb

Catalog No: #56442

Package Size: #56442-1 50ul #56442-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CD66b Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	CD66b Antibody detects endogenous levels of total CD66b
Immunogen Description	A synthesized peptide derived from human CD66b
Other Names	CD66b; CD67; CEACAM8; CGM6; NCA95;
Accession No.	Uniprot:P31997
Uniprot	P31997
Calculated MW	75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of CD66b expression in Human plasma lysate.

Product Description

Expressed in leukocytes of chronic myeloid Leukemia patients and bone marrow.

Note: This product is for in vitro research use only