

## CD83 Rabbit mAb

Catalog No: #56444

Package Size: #56444-1 50ul #56444-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	CD83 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	CD83 Antibody detects endogenous levels of total CD83
Immunogen Description	A synthesized peptide derived from human CD83
Other Names	BL11 PEN; CD83; HB15; hCD83;
Accession No.	Uniprot:Q01151
Uniprot	Q01151
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of CD83 expression in MOLT4 cell lysate.

## Product Description

May play a significant role in antigen presentation or the cellular interactions that follow lymphocyte activation.

---

Note: This product is for in vitro research use only