

## eIF4A1 Rabbit mAb

Catalog No: #56446

Package Size: #56446-1 50ul #56446-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	eIF4A1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	eIF4A1 Antibody detects endogenous levels of total eIF4A1
Immunogen Description	A synthesized peptide derived from human eIF4A1
Other Names	DDX2; DDX2A; eIF 4A I; EIF4A; eIF4A I; eIF4A-I; EIF4A1;
Accession No.	Uniprot:P60842
Uniprot	P60842
Calculated MW	46kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.


## Application Details

WB: 1:500~1:2000

ICC/IF: 1:50~1:200


FC: 1:50

## Images



Western blot analysis of eIF4A1 expression in MCF7 cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



## Product Description

ATP-dependent RNA helicase which is a subunit of the eIF4F complex involved in cap recognition and is required for mRNA binding to ribosome.

Note: This product is for in vitro research use only