

SLC6A1 / GAT1 Rabbit mAb

Catalog No: #56447

Package Size: #56447-1 50ul #56447-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SLC6A1 / GAT1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	SLC6A1 / GAT1 Antibody detects endogenous levels of total SLC6A1 / GAT1
Immunogen Description	A synthesized peptide derived from human SLC6A1 / GAT1
Other Names	GABATHG; GABATR; GABT1; GAT1; Slc6a1;
Accession No.	Uniprot:P30531
Uniprot	P30531
Calculated MW	67kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of GABA Transporter 1 / GAT 1 expression in Human fetal brain lysate.

Product Description

Terminates the action of GABA by its high affinity sodium-dependent reuptake into presynaptic terminals.

Note: This product is for in vitro research use only