

Retinoic Acid Receptor gamma Rabbit mAb

Catalog No: #56452

Package Size: #56452-1 50ul #56452-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Retinoic Acid Receptor gamma Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Rat
Specificity	Retinoic Acid Receptor gamma Antibody detects endogenous levels of total Retinoic Acid Receptor gamma
Immunogen Description	A synthesized peptide derived from human Retinoic Acid Receptor gamma
Other Names	NR1B3; RAR gamma; RARC; Rarg;
Accession No.	Uniprot:P13631
Uniprot	P13631
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 ICC/IF: 1:50~1:200

Images



Western blot analysis of Retinoic Acid Receptor gamma expression in HeLa cell lysate.

Product Description

Receptor for retinoic acid. Retinoic acid receptors bind as heterodimers to their target response elements in response to their ligands, all-trans or 9-cis retinoic acid, and regulate gene expression in various biological processes.

Note: This product is for in vitro research use only