

IRS2 Rabbit mAb

Catalog No: #56455

Package Size: #56455-1 50ul #56455-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IRS2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	IRS2 Antibody detects endogenous levels of total IRS2
Immunogen Description	A synthesized peptide derived from human IRS2
Other Names	Insulin receptor substrate 2; IRS 2; IRS-2;
Accession No.	Uniprot:Q9Y4H2
Uniprot	Q9Y4H2
Calculated MW	185kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:100

Images



Western blot analysis of IRS2 expression in HEK293 cell treated with insulin.

Product Description

May mediate the control of various cellular processes by insulin.

Note: This product is for in vitro research use only