

## Pentraxin 3 Rabbit mAb

Catalog No: #56462

Package Size: #56462-1 50ul #56462-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Pentraxin 3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse
Specificity	Pentraxin 3 Antibody detects endogenous levels of total Pentraxin 3
Immunogen Description	A synthesized peptide derived from human Pentraxin 3
Other Names	Pentaxin 3; PTX3; TNFAIP5; TSG14;
Accession No.	Uniprot:P26022
Uniprot	P26022
Calculated MW	42kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000FC:1:50

## Images



Western blot analysis of Pentraxin 3 expression in HUVEC cell lysate.

## Product Description

Plays a role in the regulation of innate resistance to pathogens, inflammatory reactions, possibly clearance of self-components and female fertility.

---

Note: This product is for in vitro research use only