

APOL1 Rabbit mAb

Catalog No: #56466

Package Size: #56466-1 50ul #56466-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	APOL1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human
Specificity	APOL1 Antibody detects endogenous levels of total APOL1
Immunogen Description	A synthesized peptide derived from human APOL1
Other Names	ApoL; APOL1; APOLI; FSGS4;
Accession No.	Uniprot:O14791
Uniprot	O14791
Calculated MW	44kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200IP:1:50FC:1:50

Images



Western blot analysis of APOL1 expression in Human plasma cell lysate.

Product Description

May play a role in lipid exchange and transport throughout the body. May participate in reverse cholesterol transport from peripheral cells to the liver.

Note: This product is for in vitro research use only