

## TBLR1 Rabbit mAb

Catalog No: #56468

Package Size: #56468-1 50ul #56468-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TBLR1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	TBLR1 Antibody detects endogenous levels of total TBLR1
Immunogen Description	A synthesized peptide derived from human TBLR1
Other Names	C21; DC42; IRA1; TBL1XR1;
Accession No.	Uniprot:Q9BZK7
Uniprot	Q9BZK7
Calculated MW	56kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

## Images



Western blot analysis of TBLR1 expression in K562 cell lysate.

## Product Description

F-box-like protein involved in the recruitment of the ubiquitin/19S proteasome complex to nuclear receptor-regulated transcription units. Plays an essential role in transcription activation mediated by nuclear receptors.

---

Note: This product is for in vitro research use only