

HNF1 beta Rabbit mAb

Catalog No: #56470

Package Size: #56470-1 50ul #56470-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HNF1 beta Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	HNF1 beta Antibody detects endogenous levels of total HNF1 beta
Immunogen Description	A synthesized peptide derived from human HNF1 beta
Other Names	FJHN; HNF 1B; HNF 2; HNF1 beta; HNF1B; HNF2; HPC11; LFB3; MODY5; TCF 2; TCF2; vHNF1;
Accession No.	Uniprot:P35680
Uniprot	P35680
Calculated MW	70kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of HNF1 beta expression in SK-OV-3 cell lysate.

Product Description

Transcription factor, probably binds to the inverted palindrome 5'-GTTAATNATTAAC-3'.

Note: This product is for in vitro research use only