

UAP1 Rabbit mAb

Catalog No: #56477

Package Size: #56477-1 50ul #56477-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	UAP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	UAP1 Antibody detects endogenous levels of total UAP1
Immunogen Description	A synthesized peptide derived from human UAP1
Other Names	AGX; AGX1; AntigenX; SPAG2; Sperm associated antigen 2; uap1;
Accession No.	Uniprot:Q16222
Uniprot	Q16222
Calculated MW	59kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50FC:1:100

Images



Western blot analysis of UAP1 expression in Jurkat cell lysate.

Product Description

Converts UDP and GlcNAc-1-P into UDP-GlcNAc, and UDP and GalNAc-1-P into UDP-GalNAc. Isoform AGX1 has 2 to 3 times higher activity towards GalNAc-1-P, while isoform AGX2 has 8 times more activity towards GlcNAc-1-P.

Note: This product is for in vitro research use only