

SRP54 Rabbit mAb

Catalog No: #56479

Package Size: #56479-1 50ul #56479-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SRP54 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP FC
Species Reactivity	Mouse, Rat, Human
Specificity	SRP54 Antibody detects endogenous levels of total SRP54
Immunogen Description	A synthesized peptide derived from human SRP54
Other Names	SRP54;
Accession No.	Uniprot:P61011
Uniprot	P61011
Calculated MW	54kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50FC:1:50

Images



Western blot analysis of SRP54 expression in 293T cell lysate.

Product Description

Binds to the signal sequence of presecretory protein when they emerge from the ribosomes and transfers them to TRAM (translocating chain-associating membrane protein).

Note: This product is for in vitro research use only