

Trefoil Factor 3 Rabbit mAb

Catalog No: #56482

Package Size: #56482-1 50ul #56482-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Trefoil Factor 3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	Trefoil Factor 3 Antibody detects endogenous levels of total Trefoil Factor 3
Immunogen Description	A synthesized peptide derived from human Trefoil Factor 3
Other Names	hITF; ITF; mITF; P1B; TFF3; TFI; TREFOIL;
Accession No.	Uniprot:Q07654
Uniprot	Q07654
Calculated MW	10kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:100

Images



Western blot analysis of Trefoil Factor 3 expression in Human colon cancer lysate.

Product Description

Involved in the maintenance and repair of the intestinal mucosa. Promotes the mobility of epithelial cells in healing processes (motogen).

Note: This product is for in vitro research use only